



- **Prodej a servis analytické  
přístrojové techniky**
- **Komplexní řešení problémů  
prvkové analýzy**

**Pernštýnská 116  
533 41 Lázně Bohdaneč  
Tel.: 466921885  
Fax.: 466921404  
E-mail: [sale@rmi.cz](mailto:sale@rmi.cz)  
Web: [www.rmi.cz](http://www.rmi.cz)**

## Nabídka firmy RMI v oblasti LIBS a LA

- [www.rmi.cz](http://www.rmi.cz)
- Kompletní systémy LIBS (LEA 500, PORTA LIBS 2000)
- Spektrometrické systémy pro LIBS  
kompletní sortiment od „low cost“  
systémů po „Hi Res“ Echelle systémy
- Lasery (Dvoupulzní ns Nd:YAG lasery)
- Detektory (intenzifikované CCD kamery)

# LEA 500

The next generation of


# LIBS

# Laser Elemental Analyzer

## LEA-S500



# Základní optické schéma systému



# **Laser – Nd:YAG dvoupulzní laser**

**Vlnová délka: 1024 nm**

**Energie pulzu: nastavitelná 80 – 150 mJ**

**Prodleva mezi pulzy: nastavitelná 0 – 20 us**

**Opakovací frekvence: 20 Hz**

**Stabilita energie pulzů: lepší jak +/- 3%**

**Chlazení: autonomní**

**Polarizace: lineární**

**Polarizace následných pulzů: kolmá**

# Laser – Nd:YAG dvoupulzní laser






Fig.2a. Single-pulse mode of laser operation



# Spektrometrický systém








## Specifikace

Focal length (mm):	500	500	500
Grating (grooves/mm):	1800	2400	3600
Wavelength range (nm):	170-800	170-600	170-400
Dispersion (nm/mm):	1.0	0.7	0.5
Wavelength resolution (nm):	0.028	0.020	0.014
Ar filled optics for wavelength range to 170 nm			

# Detekční systém – nutná vysoká rychlost – nejlépe 1 us



# **Detekční systém – specifikace**

**„Back Thinned“ and „Front Illuminated CCD“**

**2048x14 pixel**

**Spektrální rozsah 180 – 1100 nm**

**Doba načítání signálu (min.) – 1 us !!!**

**ADC: ultra fast 14 bit**

**Komunikace: Ethernet 100**



The next generation of Laser Analyzers

# Vzorkový prostor



# **Vzorkový prostor, fokusační optika**

**Evakuovaný, možnost práce v Ar atmosféře**

**Velikost erozního spotu: 0.03 – 1.2 mm**

**Zobrazení povrchu videokamerou: 1.2 x 1.2 mm**

**Alignement laser: 1 mW, 650-680 nm**

**Max. velikost vzorku: 75 x 75 x 40 mm**

**Rozsah X – Y posunu: 5 x 5 mm**

**Min. krok posunu: 1  $\mu$ m**

The next generation of Laser Analyzers


# Analyzovatelné materiály



*Cu-slítiny a čistá Cu báze*



*Fe-slítiny (surové Fe, ocel),  
Al-slítiny a čistá Al báze*



*Sklo, monokrystaly, drahé  
kameny, ...*



*Keramika*



*Plasty*

## Další typické aplikace

- ✓ *stopová analýza čistých materiálů*
- ✓ *lisované prášky*
- ✓ *vodivé & nevodivé materiály*
- ✓ *analýza vrstev, tenké vrstvy, tloušťka vrstev*

*The next generation of Laser Analyzers*

# **SOFTWARE**

***Win9x/2000/XP***



The next generation of Laser Analyzers



# Analytické Programy

Analytical programs

ProSnea - running in ADVANCED mode

Databases Tools Analyze

Scan in region

Programs Hardware Scanned spectra Calibration Samples Event log Video

**Analytical programs**

Bases: SYSTEM  
Cu  
Ti  
Fe  
SYSTEM2  
Al  
N

Programs:

- <2 FastSteel type Wx>
- <3 FastSteel type WxMoxCox>
- <1 C-steel>
- <Carbon in steel>
- <5 FeNiCo alloy Ni29Co18>
  - 51 300-330 nm Co Al Ni CuTi
  - 52 270-300 nm Si Mn
- <4 NiCrTiAl alloy typeNXTU>
- <6 Stainless steel type 12Cr18N

Select

**Active program**

Name: <1\_Glass NaCaSi>

Range: 1 - 0

**Active region**

Name: 1g261-291nm Fe Cr Mg Al

Calibration calc. method: Absolute

**Peak detection params:**

Min.Level: 0,00000 Sensitivity: 0

Store current

Sample: Glass \ BAM-5004 Program: <1\_Glass NaCaSi> Region: 1g261-291nm Fe Cr Mg Al

## The next generation of Laser Analyzers

## Control of System Status



## Základní nastavení

**ProSpectra - running in ADVANCED mode**

Databases Tools Analyze

Scan in region [v] [Play] [T] [Gears] [Color Wheel] [Copy] [Paste] [Print] [Calculator] [Help]

**Programs | Hardware | Scanned spectra | Calibration | Samples | Event log | Video**

**Region parameters | Current device state | Global parameters**

---

### Detector

Exposure time: [10] msec  
Time between scan: [10] msec  
Array for scanning: [6] [v] ☐ As active graph

**Cycle**  
☒ NonStop    Counter limit: [1] [v]  
☒ Counter  
☐ Store Series

**BackGround**  
☒ Auto Substraction  
[Get]

☐ Update View

### Laser

**Lamp parameters**  
Mode: [Internal] [v]  
Frequency: [19] [v] Hz  
Lamp energy: [5] [v] J  
Duty arc time: [30] [v] sec

**Q-SW parameters**  
Mode: [Internal] [v]  
*First channel*  
Q-SW delay: [147] [v] μs  
*Second channel*  
Q-SW delay: [153] [v] μs  
☐ Single channel

[OFF] [Q-SW] [Shutter]  
[Lamp] [Q-SW2] [Shutter2]

---

### Spectrometer

**Shutter**  
[Open] [Close]

**Filter**  
[v]

**Wavelength correction**  
WL 780.000nm L:28 C:2 R:0  
☒ use correction  
☒ auto [Change...]

**Slit width**  
[35] μm

**WaveLength**  
Current: [770] nm  
Dispersion: [0.74] nm/mm

**Turn Slit**  
[13] degree  
☒ Auto Control

### Vacuum Pump

☐ Auto Vacuum Time, msec: [7000]  
☒ Auto Clean  
☒ MaxLeak is Closed Current: [0]  
[PUMP ON] [PUMP OFF]

**Leak**  
Value: [100] [v] [Set]


[Query] [Set] [Set this params to active region]

Region: 2 280-310 nm Al Fe Mn Si Ni Sn



# Záznam spekter

Analytical programs





# Kalibrace

Analytical programs

ProSpect - running in ADVANCED mode

Databases Tools Analyze

Scan in region

Programs Hardware Scanned spectra Calibration Samples Event log Video

**Global parameters:**

☒ Calibration params  
☐ Recalibration params

Gr.Deviation, %: 2,5921  
Delta C avg, %: 1,9301  
Delta I, %: 0

**Channels:**

Al - 1  
Fe - 1  
**Sn - 1**  
Mn - 1  
Si - 1  
Ni - 1  
Al - 2

Approximation type: Line  
Analytic Line: Sn 303,4120

**Calibration Recalibration**

Standards:

[2/2] Cu alloy set 49 / 491  
[2/2] Cu alloy set 49 / 492  
[2/2] Cu alloy set 49 / 493  
[2/2] Cu alloy set 49 / 494  
[0/2] Cu alloy set 49 / 495

30.12.1899 00:00:00  
0

0,506

N 1 X 0.07 Perce 0.19

* Standard	Ct, %	CC, %	dC, %	Intensity	dIntensity	RMSD
494	0,0700	0,0680	2,908	0,1871	0	14,3997
494-1	0,0700	0,0749	10,111	0,1934	3,3926	
494-2	0,0700	0,0611	-10,176	0,1807	-3,3969	

Sample: Glass \ BAM-5004 Program: <AlFeBronze set 49> Region: 2 280-310 nm Al Fe Mn Si Ni Sn

The next generation of Laser Analyzers



Samples  
Database

# Experimental Results of Observation Field

ProSnea - running in ADVANCED mode

Databases Tools Analyze

Scan in region

Programs Hardware Scanned spectra Calibration Samples Event log Video

Sample management:

Sample description:

Groups:

Aluminium

Names:

- ☐ testarea
- ☐ 7
- ☐ 4
- ☐ 21
- ☐ 211
- ☐ 212
- ☐ 34
- ☐ 41
- ☐ 42
- ☐ 43
- ☐ 44
- ☒ 26
- ☐ test
- ☐ 126

Point analysis Area analysis

Dates:

- ☒ 23.03.2004

Programs:

- Alternative
  - Scan-1 (5x5) - [14:42]
  - Scan-2 (5x5) - [14:46]
  - Scan-3 (5x5) - [14:51]
  - Scan-4 (20x20) - [14:56]

Experiments results: ☐ view considering error

Search of analogs

	N	Cu-1	Cu-2	Fe-1	Ni-1	Si-1
2		0.1827	<0	<0	0.0964	>0
3		0.2952	<0	<0	0.0986	14.9329
4		0.4083	<0	<0	0.1834	>0
5		0.8948	0.1993	0.5128	0.5623	14.6045
6		2.1999	2.2843	1.0421	1.3106	10.6242
7		0.5844	1.7356	0.3814	0.4345	14.2429
8		1.9461	1.6049	0.6181	1.1855	11.2064
9		0.4717	<0	0.5372	0.4496	15.3522
10		0.3795	<0	0.3328	0.2988	>0

Sample: Aluminium 126 Program: <Алюминиевые литейные сплавы> Region: 275 Fe Mg Mn Si



The next generation of Laser Analyzers



Video

## Settings of Observation Field (XY stage, surface)

ProSnea - running in ADVANCED mode

Databases Tools Analyze

Scan in region

Programs Hardware Scanned spectra Calibration Samples Event log Video

**Video control panel**

Start Fit to window

Source Shots dir

Display Screenshot

Format Reset counter

Screenshot dir: Last screenshot: ImageXXXXX

**Light** 255 **Targeter** On OFF

**Telescope**

Value: 300  $\mu\text{m}$  << >>

Step: 1  $\mu\text{m}$  Calibration ...

**XY-Stage**


X: 3600  $\mu\text{m}$  DX: 100  $\mu\text{m}$

Y: -161  $\mu\text{m}$  DY: 1500  $\mu\text{m}$

Set Origin Set Delta

0:0

Origin: (0:0)



Sample: Glass \ BAM-5004 Program: <AlFeBronze set 49> Region: 2 280-310 nm Al Fe Mn Si Ni Sn



## Selection and Setting of Main Options

**Options**

**Loading**

Language: English

Interface mode: Advanced

**Calculations**

Left peak distance: 4

Right peak distance: 4

Calibration admittance: 30

**Periodic central pixel calibration**

☒ Do calibration on start

☒ Do calibration each: 40 min

Time to next calibration: 19 min

Max central pixel distance: 30

OK Cancel

*Method of  
Standard Additions  
(MSA)*







## Database of Standards and Analytical programs

### relative type (1/S)

Etalons database





Data type  
☒ Concentr. ☐ Thickness

Groups:  
102  
Pb-Sb alloy (620-74)  
Cast aluminium (137-143)  
Special Cast Iron  
Bor  
Aluminium Alloys (Mosco  
Phbr  
Bronse & Brass

Measure type: rel

Measure unit: %

Active etalon  
Add to routine

Database edit  
☐ Group ☒ Etalon  
    
☒ Fix etalon names  


	Name	Al	Cu	Mn	Ti	Cd	Mg	Fe	Si	Zn
1	1 val10	1	0.03705	0.00995	0.00339	0.00392	0.00091	0.00201	0.00127	0
2	2	1	0.0405	0.00442	0.0021	0.00168	0.00055	0.00034	0.00231	0
3	3	1	0.05256	0.00608	0.00224	0.00235	0.00027	0.00117	0.00155	0
4	4	1	0.05797	0.00287	0.00029	0.00086	0.00044	0.00043	0.00084	0
5	5	1	0	0	0.0013	0	0	0	0	0
6	131 d16	1	0.05113	0.0051	0.00065	0	0.01607	0.00521	0.0038	0.00282
7	132	1	0.02354	0.01006	0.00013	0	0.0092	0.00663	0.01241	0.00621
8	133	1	0.03063	0.00718	0.00139	0	0.02431	0.00353	0.00161	0.00161
9	134	1	0.04472	0.00636	0.00089	0	0.0069	0.007	0.00668	0.00345
10	135	1	0.05523	0.00237	0.00118	0	0.00398	0.0028	0.0008	0.00872
11	136	1	0.0383	0.00885	0	0	0.01953	0.00108	0.00949	0.00066
12	137	1	0.01884	0.01176	0	0	0.0305	0.00752	0.01612	0
13	181 amg	1	0	0.00031	0.00029	0	0.01418	0.00997	0.00055	0.00048
14	182	1	0.00239	0.00135	0.00156	0	0.02444	0.00093	0.00478	0.0027
15	183	1	0.00116	0.00285	0.00116	0	0.03494	0.00749	0.00654	0
16	184	1	0.00017	0.0097	0.00023	0	0.04743	0.00128	0.00106	0.0032
17	185	1	0.00151	0.00561	0.00063	0	0.06155	0.00572	0.00346	0.00082
18	186	1	0.00078	0.00839	0.00082	0	0.07337	0.00305	0.00164	0.00153
19	187	1	0.00333	0.01374	0	0	0.00625	0.00125	0.00968	0.00375
20	188	1	0.00043	0.00073	0.00187	0	0.09375	0.00165	0.00026	0
21	16 al9	1	0.00055	0.00199	0	0	0.0005	0.00542	0.09675	0.0005
22	17	1	0.00109	0.00404	0	0	0.00218	0.00808	0.0742	0.00164
23	18	1	0.00228	0.00565	0	0	0.00174	0.0114	0.06134	0.00326
24	19	1	0.0041	0.00626	0	0	0.00561	0.01544	0.04318	0.00497



# Wavelengths of Spectrum Lines

**Spectral lines database**

**Visualization**

☐ simple  
☒ full

**Sorting**

☐ unsorted  
☒ elements  
☐ priority

**Filters**

☒ WL : ( 326 ; 338 ) ☒ Ionisation: 3  
☐ RI : ( 0 ; 0 ) ☐ Priority from 0  
☒ Elements: Define...  
Mn,Fe,Ni,Si

**Editing**

New line Delete line Show on graph Add to routine

N	Elem	WL	Ionisation	Priority	Potential	Rel.Int	Description
4	Si	327.046	3	0	0	6	
5	Mn	328.749	3	0	0	100	
6	Fe	336.621	3	0	3.81	0	
7	Fe	335.657	3	0	3.81	0	



# Analytical Programs

Calibration programs

**Calibration programs management dialog**

Type: ☒ Concentration ☐ Thickness

Programm: <Steel\_big\_spot\_all>

Base: Fe

Regions:

- 397
- 344
- 332
- 290
- 252



ITM\_Cu  
Klapan  
<Concentration regions params>  
<Base regions params>  
<Steel\_big\_spot\_all>

Lines elements Calibration


Names: **Current element: Cr - 1**

Cr - 1  
V - 1  
Ti - 1  
Cu - 1  
Ti - 2  
V - 2


Analytic Line: Cr 329.543

Compare lines:  

Fe 328.675

Macros:  AC\_2div2.MFP

Approximation type: **Square**





## Setting of Lines and Device Calibration

Calibration programs management dialog

Type: ☒ Concentration ☐ Thickness

Programm: <Steel\_big\_spot\_all>

Base: Fe Range: 19900 - 20100

Regions:

397	[icon]
344	[icon]
332	[icon]
290	[icon]
252	[icon]

Lines elements Calibration

Type: Control

Cr 329.5  
Nb 351

Control  
Repers  
Analytic  
COMPARE

Line: [icon] [icon] [icon] fixed

Basis: [icon] [icon] [icon]

Peaks params:  
Min.Level: 0  
Sensitivity: 0  
[icon] Store current

Calibration programs management dialog

Type: ☒ Concentration ☐ Thickness

Programm: <Steel\_big\_spot\_all>

Base: Fe Range: 19900 - 20100

Regions:

397	[icon]
344	[icon]
332	[icon]
290	[icon]
252	[icon]

Lines elements Calibration

Experiments: [icon] [icon] [icon] [icon] [icon]

Etalons:

[1/1] Cr\_St\_14\_36 / 141-a  
[1/1] Cr\_St\_14\_36 / 142-a  
[1/1] Cr\_St\_14\_36 / 143-a  
[1/1] Cr\_St\_14\_36 / 144  
[1/1] YГ 0Б - YГ 9Б: YГ 0б - YГ 9б / YГ 0б  
[0/0] YГ 0Б - YГ 9Б: YГ 0б - YГ 9б / YГ 2б  
[0/1] YГ 0Б - YГ 9Б: YГ 0б - YГ 9б / YГ 3б  
[0/1] YГ 0Б - YГ 9Б: YГ 0б - YГ 9б / YГ 4б  
[1/1] YГ 0Б - YГ 9Б: YГ 0б - YГ 9б / YГ 5б  
[1/1] YГ 0Б - YГ 9Б: YГ 0б - YГ 9б / YГ 6б  
[1/1] YГ 0Б - YГ 9Б: YГ 0б - YГ 9б / YГ 7б  
[0/1] YГ 0Б - YГ 9Б: YГ 0б - YГ 9б / YГ 8б  
[0/1] YГ 0Б - YГ 9Б: YГ 0б - YГ 9б / YГ 9б



# Selection of Approximation and Calibration Curve Type

**Calibration programs management dialog**

Type: ☒ Concentration ☐ Thickness

Programm: <Steel\_big\_spot\_all>

Base: Fe Range: 19900 - 20100

Regions:

397	
344	
332	
290	
252	


Lines elements Calibration

Names: Current element: Cu - 1

Cr - 1	Analytic Line: Cu 327.396
V - 1	Compare lines:
Ti - 1	
Cu - 1	Fe 326.5617
Ti - 2	
V - 2	

Macros: AC\_2div2.MFP

Approximation type: **Line**








# Calibration Curve for C in Steel

Calibration programs





# Samples Database

Samples Database

**Samples database**

**Sample management:**

Icons: [Folder], [Folder], [Delete], [Delete], [Print], [Print], [Copy], [Copy]

**Groups:**

Common

Common  
ITM\_Cu  
Test  
Test1  
115-4

☐ H34153-1  
☐ H34153-2  
☐ H34069-2  
☐ H36442-2  
☐ 933  
☐ Andrew1  
☐ Andrew2  
☐ Andrew3  
☐ 34215m  
☐ H 34154-1  
☐ H 34154 - 2

**Sample description:**

**One point scans** | **Area scans**

**Dates:** ☐ 19.09.2002


**Programs:** ☐ <Steel\_big\_spot\_all> (7)

**Experiments results:** ☐ view considering error

N	Time	Cr-1	W-1	Al-1	Ni-1	Ni-2	Mn-1	Cr-2	V-1	Ti-1
Average	-	13.9438	0.0289	0.0569	0.2252	0.2601	0.6551	10.6176	13.1575	0.0138
Error, %		0	0	0	0	0	0	1.1025	1.3493	0.0036
1	14:21	13.9438	0.0289	0.0569						
2	14:22				0.2252	0.2601	0.6551			
3	14:22							10.1557	12.7782	0.0126
4	14:24							10.4012	12.6902	0.0128

*The next generation of Laser Analyzers*


## **Sample of Aluminum Cast Alloy**








## Visualization of sample surface, selection of surface zone to analyze




*Surface zone  
of aluminum  
cast alloy*




The next generation of Laser Analyzers

## Spectrum detection and analysis




## The next generation of Laser Analyzers

### Concentration calculation via calibration curves



## Graphic presentation of results



*element distribution  
inhomogeneity in  
aluminum cast alloy*

# PORTA LIBS 2000




EPP2000 LSR Model	Wavelength Range in nm	Grating g/mm	Slit-14 nm res.
UV2	200-400	2400	0.20
VIS4	400-600	2400	0.20
VIS4b	600-800	2400	0.20
NIR3b	800-1000	1800	0.20


EPP2000 LHR model	Wavelength Range in nm	Grating g/mm	Slit-7 nm res.
UV3	200-300	1800	0.10
UV3b	300-400	1800	0.10
VIS3	400-500	1800	0.10
VIS3b	500-600	1800	0.10
VIS3c	600-700	1800	0.10
NIR3	700-800	1800	0.10
NIR2	800-925	1200	0.10



# Solar TII optické systémy



# Solar TII optické systémy



# PROSCAN – detekční kamery





# LOTIS TII Lasery

**LOTIS TII**

/// Nd:YAG Lasers

## LS-2130 Pulsed Nd:YAG Laser



# LOTIS TII *Dvoupulzní lasery*

**LOTIS TII**

/// Nd:YAG Lasers

## LS-2134D Double Pulsed Nd:YAG Laser



# LOTIS TII *Dvoupulzní dvoubarevné*

## LS-2134D-C4 Double Pulsed Nd:YAG Laser



### Specification

Parameter		Value
Energy, mJ	1064 / 532 / 266 nm	200 / 100 / 30
Pulse duration (FWHM, at 1064 nm), ns		14–16
Pulse repetition rate, Hz		1–10
Beam divergence, mrad		≤1
Beam diameter, mm		6,3
Delay between laser pulses*, μs		1–80
Jitter**, ns		±1.0
Energy stability*** (1064 nm), %		±3.0
Size L x W x H, mm (Weight, kg)	Laser head	755 x 270 x 113 (21.0)
	Power supply	363 x 364 x 192 (16.5)
	Cooling system	363 x 364 x 192 (15.5)
	Remote control	105 x 175 (0.5)
Power requirements		Single phase, 220±20 V, 50–60 Hz, 600 W

\* 1 μs steps, other steps (1–100 ns) are available on request

\*\* with respect to external trigger of Q-switch

\*\*\* shot to shot for 99% of pulses

# LOTIS TII *Dvoupulzní lasery*

## *s vyšší opakovací frekvencí – 50 Hz*

**LOTIS TII**

/// Nd:YAG Lasers

**LS-2136D-5** Double Pulsed  
High Repetition Rate Nd:YAG Laser



### Specification

Parameter	Value	
Energy, mJ	1064 / 532 nm	105 / 52
Pulse duration (FWHM, at 1064 nm), ns	14–16	
Pulse repetition rate, Hz	1–50	
Beam divergence, mrad	≤1	
Beam diameter, mm	5	
Delay between laser pulses*, μs	1–80	
Jitter**, ns	±1.0	
Energy stability*** (1064 nm), %	±3.0	
Size L x W x H, mm (Weight, kg)	Laser head	755 x 270 x 113 (23.0)
	Power supply	446 x 449 x 177 (20.0)
	Cooling system	446 x 449 x 266 (23.0)
	Remote control	105 x 175 (0.5)
Power requirements	Single phase, 220±20 V, 50–60 Hz, 2000 W	

\* 1 μs steps, other steps (1–100 ns) are available on request

\*\* with respect to external trigger of Q-switch

\*\*\* shot to shot for 99% of pulses

# LOTIS TII Lasery

**LOTIS TII**

/// Nd:YAG Lasers

## LS-2134UTF Pulsed Nd:YAG Laser with VRM and Built-in Harmonics



### Specification

Parameter		Value
Energy, mJ	1064 / 532 / 355 / 266 nm	270 / 170 / 60 / 40
Pulse duration, ns	1064 / 532 / 355 / 266 nm	7–8 / 6–7 / 5–7 / 5–6
Pulse repetition rate, Hz		1–15
Beam divergence, mrad		≤0.8
Beam diameter, mm		≤6.0
Jitter*, ns		±1.0
Pointing stability, mrad		0.1
Energy stability** (1064 nm), %		±3.0
Size L x W x H, mm (Weight, kg)	Laser head	653 x 185 x 113 (14.0)
	Power supply	363 x 364 x 192 (15.5)
	Cooling system	363 x 364 x 192 (12.5)
	Remote control	105 x 175 (0.5)
Power requirements		Single phase, 220±20 V, 50–60 Hz, 750 W

\* with respect to external trigger of Q-switch

\*\* shot to shot for 99% of pulses

# LOTIS TII Lasery 213 nm

**LOTIS TII**

/// Nd:YAG Lasers

**LS-2134UF-5** Pulsed Nd:YAG Laser with VRM and Built-in Harmonics



## Specification

Parameter		Value
Energy, mJ	532 / 213 nm	170 / 6
Pulse duration, ns	532 / 213 nm	6–7 / 5–6
Pulse repetition rate, Hz		1–15
Beam divergence, mrad		≤0.8
Beam diameter, mm		≤6.0
Jitter*, ns		±1.0
Pointing stability, mrad		0.1
Energy stability** (1064 nm), %		±3.0
Size L x W x H, mm (Weight, kg)	Laser head	653 x 185 x 113 (14.0)
	Power supply	363 x 364 x 192 (15.5)
	Cooling system	363 x 364 x 192 (12.5)
	Remote control	105 x 175 (0.5)
Power requirements		Single phase, 220±20 V, 50–60 Hz, 750 W


\* with respect to external trigger of Q-switch

\*\* shot to shot for 99% of pulses



# Solar TII Lasery

## WIDELY-TUNABLE Ti:SAPPHIRE LASER




Widely-Tunable  
Ti:Sapphire Laser

MODELS

**CF 131A**  
**CF 131MA**

- Widely-tunable from IR to UV (950-189 nm)
- All solid-state
- Narrow and super narrow linewidth
- Pollution-free
- PC control (Model CF131A)



# LOTIS TII - lasery

## WIDELY-TUNABLE PULSED Ti:SAPPHIRE LASER



Widely-Tunable  
Pulsed Ti:Sapphire Laser

MODEL  
CF 125

- High efficiency
- Built-in SHG
- Wide tunable range with single cavity optics set

### TUNABLE LASER PARAMETERS

Active medium:	Ti:Sapphire
Tuning range	
at fundamental:	690-1000 nm
at second harmonic:	350-500 nm
Pump conversion efficiency in the maximum of tuning curve	
at fundamental:	$\leq 25\%$
at second harmonic:	$\leq 8\%$
Linewidth:	0,1-0,2 nm
Beam divergence:	$< 1,5$ mrad
Overall size (LxWxH):	538x150x72 mm
Weight:	6,5 kg

### PUMP PARAMETERS

Wavelength:	532 nm
Pulse energy (max):	300 mJ
Pulse duration:	8-20 ns
Power (max):	15 W
Pulse repetition rate (max):	50 Hz
Beam diameter (max):	9 mm

# LOTIS TII - lasery

## QUASI-CW Ti:SAPPHIRE LASER



Quasi-CW  
Ti:Sapphire Laser

MODEL  
**CF 231**

- KHz pulse repetition rate
- Narrow linewidth
- Wide spectral range
- Built-in SHG

- ◆ High-quality output is generated in a computer-optimized aberration-free cavity. Low threshold generation is provided by extra-quality optics in use.
- ◆ The output coupler is adapted to variation of the Ti:Sapphire gain. A single set of cavity optics provides the extremely wide 695-960 nm tuning range.
- ◆ The built-in SHG unit with cylindric optics and an efficient nonlinear crystal extend the tunability range from 350 to 475 nm.

# LOTIS TII - lasery

## SOLID STATE RAMAN SHIFTER



Solid State  
Raman Shifter

MODEL

**CRS 14**

**CRS14 SOLID STATE RAMAN SHIFTER FOR CONVERSION OF  
Nd:YAG LASER OUTPUT**

# LOTIS TII - lasery

## AUTOMATIC OPTICAL PARAMETRIC OSCILLATOR



Extremely Widely-Tunable  
Automatic Optical  
Parametric Oscillator

MODEL  
**COPO 2200A**

- 210-340 nm/420-690 nm/730-2200 nm tuning range
- Automatic operation and wavelength control
- Availability of adaptation to a variety of pump lasers

# LA-0TOF-MS




# OptiMass 9500






## Kompletní spektrum rozloženého vzorku skla.






## Kompletní spektrum rozloženého vzorku skla.





# Co je v tomto vzorku ?





# Prvkový „otisk prstu“

